

PRESSALIT.

Product	Waste basket, 297 x 222 x 344 mm
Item number	RT764
	VVS:615801764
	NRF:7501031
	EAN:5708590207858
	HMI: 124561
	NBS:N13:449
Product description	Waste basket, A:297 mm X B:222 mm x C:344 mm.
	Waste basket, A:297 mm X B:222 mm x C:344 mm.
Materials	Basic material: Steel
	Surface material: Rilsan
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for



PRESSALIT.®

Product	Waste basket, 297 x 222 x 344 mm
Item number	RT764
Cleaning	Baskets. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use corrosive or abrasive cleaners as it may result in the surface being damaged. For surface disinfection, chlorine-based disinfectants may be used according to the product specification and safety precautions for these products. The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container. DISINFECTANTS: Clinimax Difficil-S Erisan OXY+ Sprit (husholdningssprit) Suma Bac D10 Virkon Rodalon Desinfektionsspray Lysoformin Spezial CaviCide Sufanios Premium Tvättsprit 70% WipeClean UniversalCleaning Chrolorine Disinfection Wipes Klorilli (Kiilto Clean) NON-DISINFECTANTS: Suma Star free D1 Taski 100 Taski 200 Milizid Citro The product must be wiped down with a soft rag wrung out in clear water immediate after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.



PRESSALIT®

Product Waste basket, 297 x 222 x 344 mm

Item number RT764

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com