




PRESSALIT®

Product Item number	Diaper basket, 402 x 212 x 174 mm RT767
Product description	<p>VVS:615801767 NRF:7501034 EAN:5708590171074 HMI: 124570 NBS:N13:449</p>  <p>Basket for nappies, A:402 mm x B:212 mm x C: 174 mm. Basket for nappies, A:402 mm x B:212 mm x C: 174 mm.</p>
Materials	<p>Basic material: Steel</p> <p>Surface material: Rilsan</p>
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



PRESSALIT®

Product	Diaper basket, 402 x 212 x 174 mm
Item number	RT767
Cleaning	<p>Baskets.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>For surface disinfection, chlorine-based disinfectants may be used according to the product specification and safety precautions for these products.</p> <p>The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container.</p> <p>DISINFECTANTS:</p> <p>Clinimax Difficil-S Erisan OXY+ Sprit (husholdningsprit) Suma Bac D10 Virkon Rodalon Desinfektionsspray Lysoformin Spezial CaviCide Sufanios Premium Tvättsprit 70% WipeClean UniversalCleaning Chlorine Disinfection Wipes Klorilli (Kiilto Clean)</p> <p>NON-DISINFECTANTS:</p> <p>Suma Star free D1 Taski 100 Taski 200 Milizid Citro</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p>



PRESSALIT®

Product	Diaper basket, 402 x 212 x 174 mm
Item number	RT767

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®